Dragon Medical One – PowerMic Mobile App

SETTING UP POWERMIC MOBILE APP

Follow the steps below to set up the PowerMic Mobile app to allow you to use your smart phone as a wireless PowerMic with Dragon Medical One.

- 1. Download the **PowerMic Mobile** app on your iPhone or Android but you don't need to open it. If prompted on Android to Review Google Account, choose SKIP. If prompted to allow access to the microphone click YES.
- 2. You need to access your **Welcome Kit PDF** on your phone.
- 3. Open the Welcome Kit on your phone.

Scroll to the page: MOBILE APP INSTALLATION & CONFIGURATION.

4. Under STEP2 on that page:

Click the link that shows either "Click to configure iOS" for iPhone users or "Click here to configure Android", depending on your phone.

- 5. If prompted, choose "Open with the PowerMic Mobile app" or it may automatically open in the PowerMic Mobile app. Accept Agreement and click "Got it"
- 6. Enter a profile name (or accept default of Profile1 or 2) and then enter your Dragon Medical One user (your email address that you log into DMO with) and login.

TROUBLESHOOTING

if you are having issues with the configuration link and it will not open PowerMic Mobile then you can connect manually by pasting the URL directly into the App.

You can paste the URL link from your Welcome Kit into your phone where it shows Paste your organization profile URL here.

Here is a copy of each URL link below:

iPhone URL:

dmic://config_?NmsToken=NjhBMUI3MDItN0RCNy00QzYzLTlERTktOEQwNTJDNzZDOUUw&NmsBaseUrl=nms-ca.nuancehdp.com

Android URL:

http://config_?NmsToken=NjhBMUI3MDItN0RCNy00QzYzLTIERTktOEQwNTJDNzZDOUUw&NmsBaseUrl=nms-ca.nuancehdp.com

Note: To use the PowerMic mobile microphone **FIRST log into Dragon Medical One** with your username and make sure Microphone is set to PowerMic Mobile and **THEN launch the PowerMic Mobile** app otherwise it will show "searching for computer".

If you require assistance with installing Dragon Medical One, you are welcome to contact Ice Peak Solutions.